

Product Datasheet

Dog IL1RA ELISA Kit (orb403466)

Description	Dog Interleukin-1 receptor antagonist protein(IL1RN) ELISA kit
Reactivity	Canine
Range	15.6 pg/ml-1000 pg/ml
Target	IL1RN
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9BEH0
Sensitivity	3.9 pg/ml
Expiration Date	Please enquire.