

Product Datasheet

Dog S100A12 ELISA Kit (orb403508)

Description	Canine S100 calcium binding protein A12/Calgranulin-C(S100A12) ELISA Kit
Reactivity	Canine
Range	0.625 ng/mL-40 ng/mL
Target	S100A12
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	F1PMK3
Sensitivity	0.156 ng/mL
Expiration Date	Please enquire.