

Product Datasheet

Dog IL1B ELISA Kit (orb403520)

Description	Dog Interleukin 1 β (IL- 1 β) ELISA Kit
Reactivity	Canine
Range	7.8 pg/mL-500 pg/mL
Target	IL1B
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q28292
Sensitivity	1.95 pg/mL
Expiration Date	Please enquire.