

## Product Datasheet

### Dog FSH ELISA Kit (orb403535)

<b>Description</b>	Canine follicle-stimulating hormone,FSH ELISA Kit
<b>Reactivity</b>	Canine
<b>Range</b>	0.35 mIU/mL-140 mIU/mL
<b>Target</b>	FSH
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<b>E2RJZ8</b>
<b>Sensitivity</b>	0.35 mIU/mL
<b>Expiration Date</b>	Please enquire.