

Product Datasheet

Dog MMP-13 ELISA Kit (orb403553)

Description	Canine Matrix metalloproteinase 13,MMP-13 ELISA kit
Reactivity	Canine
Range	0.156 ng/mL-10 ng/mL
Target	MMP-13
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	F1PVS9
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.