

Product Datasheet

Dog MMP2/Gelatinase A ELISA Kit (orb403554)

Description	Canine Matrix metalloproteinase 2/Gelatinase A,MMP-2 ELISA kit
Reactivity	Canine
Range	4.69 ng/mL-300 ng/mL
Target	MMP2/Gelatinase A
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	F1PMK7
Sensitivity	1.17 ng/mL
Expiration Date	Please enquire.