

Product Datasheet

Dog MMP3 ELISA Kit (orb403555)

Description Dog Matrix metalloproteinase 3,MMP-3
ELISA kit

Reactivity Canine

Range 1.56 ng/mL-100 ng/mL

Target MMP3

Note For research use only

Application notes Detection Wavelength: 450 nm

Sample Types serum, plasma, cell culture supernates,
urine, tissue homogenates

Uniprot ID **Q6Y4Q5**

Sensitivity 0.39 ng/mL

Expiration Date Please enquire.