

## Product Datasheet

### Dog MMP9 ELISA Kit (orb403556)

<b>Description</b>	Dog matrix metalloproteinase 9/Gelatinase B, MMP-9 ELISA Kit
<b>Reactivity</b>	Canine
<b>Range</b>	12.5 ng/mL-800 ng/mL
<b>Target</b>	MMP9
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma
<b>Uniprot ID</b>	<b>O18733</b>
<b>Sensitivity</b>	3.12 ng/mL
<b>Expiration Date</b>	Please enquire.