

## Product Datasheet

### Dog RV antibody (IgG) ELISA Kit (orb403562)

<b>Description</b>	Dog rabies virus (RV) antibody (IgG) ELISA Kit
<b>Reactivity</b>	Canine
<b>Range</b>	Request Information
<b>Target</b>	RV Ab (IgG)
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum
<b>Sensitivity</b>	Request Information
<b>Expiration Date</b>	Please enquire.