

## Product Datasheet

### Fish IL8/CXCL8 ELISA Kit (orb403637)

|                          |                                   |
|--------------------------|-----------------------------------|
| <b>Description</b>       | Fish Interleukin-8,IL-8 ELISA Kit |
| <b>Reactivity</b>        | Fish                              |
| <b>Range</b>             | 31.25 pg/mL-500 pg/mL             |
| <b>Target</b>            | IL8/CXCL8                         |
| <b>Note</b>              | For research use only             |
| <b>Application notes</b> | Detection Wavelength: 450 nm      |
| <b>Sample Types</b>      | serum, plasma                     |
| <b>Uniprot ID</b>        | <b>Q6F4N2</b>                     |
| <b>Sensitivity</b>       | 31.25 pg/mL                       |
| <b>Expiration Date</b>   | Please enquire.                   |