

Product Datasheet

Fish PRL/LTH ELISA Kit (orb403641)

Description	Fish Prolactin/Luteotropic Hormone(PRL/LTH) ELISA Kit
Reactivity	Fish
Range	10 µIU/mL-4000 µIU/mL
Target	PRL/LTH
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P09319
Sensitivity	10 µIU/mL
Expiration Date	Please enquire.