

Product Datasheet

Fish GnRH ELISA Kit (orb403649)

Description	Fish gonadotropin-releasing hormone, GnRH ELISA Kit
Reactivity	Fish
Range	0.5 ng/mL-20 ng/mL
Target	GnRH
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates, cell lysates
Uniprot ID	Q5Y835
Sensitivity	0.5 ng/mL
Expiration Date	Please enquire.