

Product Datasheet

Goat CGA ELISA Kit (orb403695)

| | |
|--------------------------|---|
| Description | Goat Thyroid Stimulating Hormone, TSH ELISA Kit |
| Reactivity | Goat |
| Range | 0.5 µU/mL-20 µU/mL |
| Target | CGA |
| Note | For research use only |
| Application notes | Detection Wavelength: 450 nm |
| Sample Types | serum, plasma, tissue homogenates |
| Uniprot ID | Q8WMW8 |
| Sensitivity | 0.25 µU/mL |
| Expiration Date | Please enquire. |