

Product Datasheet

Goat HSPA2 ELISA Kit (orb403719)

Description	Goat Heat Shock Protein 70(HSP-70) ELISA KIT
Reactivity	Goat
Range	23.44 ng/mL-1500 ng/mL
Target	HSPA2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9TUG3
Sensitivity	5.86 ng/mL
Expiration Date	Please enquire.