

## Product Datasheet

### Goat GH1 ELISA Kit (orb403720)

|                          |                                   |
|--------------------------|-----------------------------------|
| <b>Description</b>       | Goat growth hormone,GH ELISA KIT  |
| <b>Reactivity</b>        | Goat                              |
| <b>Range</b>             | 6.25 ng/mL-100 ng/mL              |
| <b>Target</b>            | GH1                               |
| <b>Note</b>              | For research use only             |
| <b>Application notes</b> | Detection Wavelength: 450 nm      |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates |
| <b>Uniprot ID</b>        | <b>P67931</b>                     |
| <b>Sensitivity</b>       | 6.25 ng/mL                        |
| <b>Expiration Date</b>   | Please enquire.                   |