

Product Datasheet

GGCT antibody (orb40375)

Description Peptide affinity purified rabbit polyclonal antibody.

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications WB

Immunogen Synthetic Peptide

Target GGCT

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal

MW 21008

Uniprot ID [O75223](#)

NCBI [NP_001186745.1](#), [NP_076956.1](#), [NP_001186746.1](#), [NP_001186744.1](#)

Dilution Range WB: 1:1000

Expiration Date 12 months from date of receipt.

