

Product Datasheet

Guinea Pig IL10 ELISA Kit (orb403750)

Description	Guinea pig interleukin 10,IL-10 ELISA Kit
Reactivity	Guinea pig
Range	31.25 pg/mL-2000 pg/mL
Target	IL10
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9Z1Y5
Sensitivity	7.81 pg/mL
Expiration Date	Please enquire.