

Product Datasheet

Guinea Pig IL1B ELISA Kit (orb403755)

Description	Guinea pig Interleukin-1 beta(IL1B) ELISA kit
Reactivity	Guinea pig
Range	15.6 pg/mL-1000 pg/mL
Target	IL1B
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9WVG1
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.