

Product Datasheet

Guinea Pig TNF ELISA Kit (orb403782)

Description	Guinea pig Tumor necrosis factor α ,TNF- α ELISA Kit
Reactivity	Guinea pig
Range	3.12 pg/mL-200 pg/mL
Target	TNF
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P51435
Sensitivity	0.78 pg/mL
Expiration Date	Please enquire.