

Product Datasheet

Hamster IL-1 β ELISA Kit (orb403789)

Description	Hamster Interleukin 1 β , IL-1 β ELISA Kit
Reactivity	Golden Hamster
Range	12.5 pg/mL-800 pg/mL
Target	IL1B
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q9WVQ9
Sensitivity	3.12 pg/mL
Expiration Date	Please enquire.