

Product Datasheet

Horse HIF1A ELISA Kit (orb403859)

Description	Horse hypoxia inducible factor 1, alpha subunit (basic helix-loop-helix transcription factor) (HIF1A...
Reactivity	Equine
Range	31.25 pg/mL-2000 pg/mL
Target	HIF1A
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	F6VG41
Sensitivity	7.81 pg/mL
Expiration Date	Please enquire.