

Product Datasheet

Human SFN ELISA Kit (orb403874)

| | |
|--------------------------|---|
| Description | Human 14-3-3 protein sigma(SFN) ELISA kit |
| Reactivity | Human |
| Range | 0.27 ng/ml- 200ng/ml |
| Target | SFN |
| Note | For research use only |
| Application notes | Detection Wavelength: 450 nm |
| Sample Types | serum, plasma, tissue homogenates |
| Uniprot ID | P31947 |
| Sensitivity | 0.27 ng/ml |
| Expiration Date | Please enquire. |