

Product Datasheet

Human HSD17B10 ELISA Kit (orb403883)

Description Human 3-hydroxyacyl-CoA dehydrogenase type-2(HSD17B10) ELISA kit

Reactivity Human

Range 15.625 pg/ml- 1000pg/ml

Target HSD17B10

Note For research use only

Application notes Detection Wavelength: 450 nm

Sample Types tissue homogenates, cell lysates.

Uniprot ID **Q99714**

Sensitivity 3.9 pg/ml

Expiration Date Please enquire.