

Product Datasheet

Human ADAMTS18 ELISA Kit (orb403898)

Description	Human A disintegrin and metalloproteinase with thrombospondin motifs 18(ADAMTS18) ELISA kit
Reactivity	Human
Range	0.016ug/ml- 10ug/ml
Target	ADAMTS18
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q8TE60
Sensitivity	0.016ug/ml
Expiration Date	Please enquire.