

Product Datasheet

Human ADAMTS2 ELISA Kit (orb403900)

Description	Human A disintegrin and metalloproteinase with thrombospondin motifs 2(ADAMTS2) ELISA kit
Reactivity	Human
Range	3.12 ng/mL-200 ng/mL
Target	ADAMTS2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	O95450
Sensitivity	0.78 ng/mL
Expiration Date	Please enquire.