

Product Datasheet

Human ADAMTS5 ELISA Kit (orb403903)

Description	Human A disintegrin and metalloproteinase with thrombospondin motifs 5(ADAMTS5) ELISA kit
Reactivity	Human
Range	15.6 ng/mL-1000 ng/mL
Target	ADAMTS5
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9UNA0
Sensitivity	3.9 ng/mL
Expiration Date	Please enquire.