

Product Datasheet

Human ARPC2 ELISA Kit (orb403916)

Description	Human Actin-related protein 2/3 complex subunit 2(ARPC2) ELISA kit
Reactivity	Human
Range	23.5 pg/mL-1500 pg/mL
Target	ARPC2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	O15144
Sensitivity	5.8 pg/mL
Expiration Date	Please enquire.