

Product Datasheet

Human ARC ELISA Kit (orb403919)

Description	Human Activity-regulated cytoskeleton-associated protein(ARC) ELISA kit
Reactivity	Human
Range	25 pg/mL-1600 pg/mL
Target	ARC
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q7LC44
Sensitivity	6.25 pg/mL
Expiration Date	Please enquire.