

Product Datasheet

Human ADAMTSL2 ELISA Kit (orb403928)

Description	Human ADAMTS-like protein 2(ADAMTSL2) ELISA kit
Reactivity	Human
Range	3.12 ng/mL-200 ng/mL
Target	ADAMTSL2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q86TH1
Sensitivity	0.78 ng/mL
Expiration Date	Please enquire.