

Product Datasheet

Human ADAMTSL4 ELISA Kit (orb403930)

Description	Human ADAMTS-like protein 4(ADAMTSL4) ELISA kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Target	ADAMTSL4
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q6UY14
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.