

## Product Datasheet

### Human CRYAB ELISA Kit (orb403955)

|                          |                                                 |
|--------------------------|-------------------------------------------------|
| <b>Description</b>       | Human Alpha-crystallin B chain(CRYAB) ELISA kit |
| <b>Reactivity</b>        | Human                                           |
| <b>Range</b>             | 12.5 pg/mL-800 pg/mL                            |
| <b>Target</b>            | CRYAB                                           |
| <b>Note</b>              | For research use only                           |
| <b>Application notes</b> | Detection Wavelength: 450 nm                    |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates, cell lysates |
| <b>Uniprot ID</b>        | <b>P02511</b>                                   |
| <b>Sensitivity</b>       | 3.12 pg/mL                                      |
| <b>Expiration Date</b>   | Please enquire.                                 |