

Product Datasheet

Human ANKRD26 ELISA Kit (orb403976)

Description	Human Ankyrin repeat domain-containing protein 26(ANKRD26) ELISA kit
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Target	ANKRD26
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q9UPS8
Sensitivity	7.81 pg/mL
Expiration Date	Please enquire.