

Product Datasheet

Human AGR3 ELISA Kit (orb403979)

Description	Human Anterior gradient protein 3 homolog(AGR3) ELISA kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Target	AGR3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q8TD06
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.