



Product Datasheet

HFE2 antibody (orb40400)

Description	Rabbit polyclonal antibody to HFE2 also known as Hemojuvelin, is a membrane-bound and soluble protei...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	HJV
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	45080
Uniprot ID	Q6ZVN8
NCBI	NP_998817.1 , NP_660320.3 , NP_973733.1 , NP_998818.1
Dilution Range	WB: 1:1000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...