

Product Datasheet

Human BCL2A1 ELISA Kit (orb404019)

Description	Human Bcl-2-related protein A1(BCL2A1) ELISA kit
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Target	BCL2A1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q16548
Sensitivity	7.81 pg/mL
Expiration Date	Please enquire.