

Product Datasheet

Human B4GALT1 ELISA Kit (orb404021)

Description	Human Beta-1,4-galactosyltransferase 1(B4GALT1) ELISA kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Target	B4GALT1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell lysates, tissue homogenates
Uniprot ID	P15291
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.