



Product Datasheet

INSC antibody (orb40404)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	INSC
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	63469
Uniprot ID	Q1MX18
NCBI	NP_001265244.1 , NP_001027024.3 , NP_001265245.1 , NP_001036001.1 , NP_001265242.1 , NP_001265243.1
Dilution Range	WB: 1:1000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of INSC protein in WiDr cell lysate.