

Product Datasheet

Human BAMBI ELISA Kit (orb404040)

Description	Human BMP and activin membrane-bound inhibitor homolog(BAMBI) ELISA kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Target	BAMBI
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q13145
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.