

## Product Datasheet

### Human BMP3 ELISA Kit (orb404044)

**Description** Human Bone morphogenetic protein 3(BMP3) ELISA kit

**Reactivity** Human

**Range** 15.6 pg/mL-1000 pg/mL

**Target** BMP3

**Note** For research use only

**Application notes** Detection Wavelength: 450 nm

**Sample Types** serum, plasma, tissue homogenates

**Uniprot ID** [P12645](#)

**Sensitivity** 3.9 pg/mL

**Expiration Date** Please enquire.