

Product Datasheet

Human BMP8A ELISA Kit (orb404046)

Description	Human Bone morphogenetic protein 8A(BMP8A) ELISA kit
Reactivity	Human
Range	37.5 pg/mL-2400 pg/mL
Target	BMP8A
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q7Z5Y6
Sensitivity	9.38 pg/mL
Expiration Date	Please enquire.