

Product Datasheet

Human CEACAM8 ELISA Kit (orb404093)

Description	Human Carcinoembryonic antigen-related cell adhesion molecule 8(CEACAM8) ELISA kit
Reactivity	Human
Range	0.625 ng/mL-40 ng/mL
Target	CEACAM8
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P31997
Sensitivity	0.156 ng/mL
Expiration Date	Please enquire.