

Product Datasheet

Human CLCF1 ELISA Kit (orb404094)

Description	Human Cardiotrophin-like cytokine factor 1 (CLCF1) ELISA kit
Reactivity	Human
Range	0.312 ng/mL-20 ng/mL
Target	CLCF1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q9UBD9
Sensitivity	0.078 ng/mL
Expiration Date	Please enquire.