

## Product Datasheet

### Human CSN1S1 ELISA Kit (orb404099)

<b>Description</b>	Human casein alpha s1 (CSN1S1) ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	37.5 ?g/mL-600 ?g/mL
<b>Target</b>	CSN1S1
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	breast milk
<b>Uniprot ID</b>	<a href="#">P47710</a>
<b>Sensitivity</b>	18.75 ?g/mL
<b>Expiration Date</b>	Please enquire.