

## Product Datasheet

### Human CSN1S1 ELISA Kit (orb404099)

|                          |   |
|--------------------------|---|
| <b>Description</b>       | Human casein alpha s1 (CSN1S1)<br>ELISA kit |
| <b>Reactivity</b>        | Human                                       |
| <b>Range</b>             | 37.5 µg/mL-600 µg/mL                        |
| <b>Target</b>            | CSN1S1                                      |
| <b>Note</b>              | For research use only                       |
| <b>Application notes</b> | Detection Wavelength: 450 nm                |
| <b>Sample Types</b>      | breast milk                                 |
| <b>Uniprot ID</b>        | <b>P47710</b>                               |
| <b>Sensitivity</b>       | 18.75 µg/mL                                 |
| <b>Expiration Date</b>   | Please enquire.                             |