



Product Datasheet

GDA antibody (orb40412)

Description	Rabbit polyclonal antibody to GDA
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthetic Peptide
Target	GDA
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	51003
Uniprot ID	Q9Y2T3
NCBI	NP_001229434.1 , NP_001229435.1 , NP_001229436.1 , NP_004284.1
Dilution Range	WB: 1:1000, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of in HL-60 cell l...