

Product Datasheet

Human CSPG4 ELISA Kit (orb404145)

Description	Human Chondroitin sulfate proteoglycan 4(CSPG4) ELISA kit
Reactivity	Human
Range	25 pg/mL-1600 pg/mL
Target	CSPG4
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q6UVK1
Sensitivity	6.25 pg/mL
Expiration Date	Please enquire.