

Product Datasheet

Human COL10A1 ELISA Kit (orb404172)

Description	Human Collagen alpha-1(X) chain (COL10A1) ELISA kit
Reactivity	Human
Range	0.312 ng/mL-20 ng/mL
Target	COL10A1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q03692
Sensitivity	0.078 ng/mL
Expiration Date	Please enquire.