

Product Datasheet

FAM40A antibody (orb40419)

Description Peptide affinity purified rabbit polyclonal antibody.

Species/Host Rabbit

Reactivity Human, Mouse

Conjugation Unconjugated

Tested Applications FC, IHC-P, WB

Immunogen Synthetic Peptide

Target STRIP1

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal



MW 95576

Uniprot ID [Q5VSL9](#)


NCBI [NP_001257697.1](#), [NP_149079.2](#)

Dilution Range WB: 1:1000, IHC-P: 1:10-50, FACS: 1:10-50

Expiration Date 12 months from date of receipt.



Immunohistochemical analysis of formalin...



Western blot analysis of FAM40A antibody...