

Product Datasheet

Human CLEC2D ELISA Kit (orb404230)

Description	Human C-type lectin domain family 2 member D(CLEC2D) ELISA kit
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Target	CLEC2D
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q9UHP7
Sensitivity	7.81 pg/mL
Expiration Date	Please enquire.