

Product Datasheet

Human CLEC9A ELISA Kit (orb404236)

Description	Human C-type lectin domain family 9 member A(CLEC9A) ELISA kit
Reactivity	Human
Range	0.312 ng/mL-20 ng/mL
Target	CLEC9A
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q6UXN8
Sensitivity	0.078 ng/mL
Expiration Date	Please enquire.